



Dear PessaryAssistant™ Customer,

Thank you for your purchase of PessaryAssistant™, which I invented as a mechanical aid for pessary placement & removal. It is CE-certified for sale in the EU, UK, Australia, & other regions.

This new & improved version is made from medical-grade Polypropylene (PP), making it lighter & more flexible. It is tapered to allow hook retrieval through small holes in ring pessary diaphragms. And the slightly larger-diameter hook enables retrieval of the stem of Gellhorn pessaries.

Enclosed you will find your PessaryAssistant™, which includes:

- 1) A “U”-shaped (Fold) side which helps hold a folded ring pessary in preparation for placement, &
- 2) A “C”-shaped (Hook) side which helps hook & remove a ring pessary or Gellhorn stem.

Also enclosed is a corresponding number of Instructions for Use (IFU) sheets.

Additional resources, including hand-demonstration videos & IFUs in multiple languages, are available on my website. All of the following URLs link to the same site. Or you may use the QR code below:



www.PessaryAssistant.com

www.PessaryAssist.com

www.PessaryHelp.com

www.PessaryHook.com

www.PessaryTool.com

I continually refine PessaryAssistant™ based on patient & clinician feedback. Your input is always appreciated & helps guide these ongoing improvements.

Thank you again for your purchase!

Warmly,

Chris Radpour MD

cradpour@galenmedical.com